

Product Information




Palitantin

Cat No: 25208 - 500 µg

General Data

Shipping:	dry ice
Formulation:	A solid
Purity:	≥95%



Product Overview

Palitantin is a fungal metabolite originally isolated from *P. palitans*.{39756} It inhibits HIV-1 integrase with IC50 values of 350 and 268 µM for 3'-processing and strand transfer activity, respectively.{39757} Palitantin has been used as a standard for dereplication of natural products.{41781}

FP/22/24

For research laboratory use only – Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

Contact Bertin Bioreagent

Parc d'Activités du Pas du Lac
10 bis avenue Ampère
78180 Montigny le Bretonneux - France
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>
<https://www.bertin-bioreagent.com>