




Butylone-d₃ (hydrochloride) (CRM)

Cat No: 28819 - 100 µg

General Data

Shipping: dry ice

Formulation: A 100 µg/ml solution in methanol



Special advice: Reference Materials RMs and Certified Reference Materials CRMs are dispatched directly to the end-user in accordance with ISO requirements. Please check regulation status before ordering; additional fees can apply.

Product Overview

Butylone-d₃ (hydrochloride) (CRM) (Item No. 28819) is a certified reference material intended for use as an internal standard for the quantification of butylone (Item Nos. 10393 16059) by GC- or LC-MS. Butylone is categorized as a cathinone. It is a metabolite of bk-DMBDB and has been detected in products sold as bath salts. Butylone-d₃ is regulated as a Schedule I compound in the United States. Butylone-d₃ (hydrochloride) (CRM) (Item No. 28819) is provided as a DEA exempt preparation. This product is intended for research and forensic applications.

FP/24/24

For research laboratory use only – Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

Contact Bertin Bioreagent

Parc d'Activités du Pas du Lac
10 bis avenue Ampère
78180 Montigny le Bretonneux - France
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>
<https://www.bertin-bioreagent.com>