

# Product Information




## MIV-150

Cat No: 9002360 - 500 µg

### General Data

Shipping:	dry ice
Formulation:	A crystalline solid
Purity:	≥98%



### Product Overview

MIV-150 is a tight-binding, allosteric inhibitor of reverse transcriptase that is active against HIV-1 and HIV-2 (both EC50s = 1 nM in vitro).{29583} It has been shown to inactivate viruses that are resistant to other antiviral drugs, including non-nucleoside reverse transcriptase inhibitors, nucleoside reverse transcriptase inhibitors, and protease inhibitors.{29583} Although MIV-150 possesses poor oral bioavailability, it demonstrates efficacy when formulated as a topical microbicide.{29585,29584}

FP/22/24

For research laboratory use only – Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

Contact Bertin Bioreagent

Parc d'Activités du Pas du Lac  
10 bis avenue Ampère  
78180 Montigny le Bretonneux - France  
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>  
<https://www.bertin-bioreagent.com>