




Dihydrocodeine (CRM)

Cat No: ISO60142 - 1 mg

General Data

Shipping: dry ice

Formulation: A 1 mg/ml solution in methanol



Special advice: Reference Materials RMs and Certified Reference Materials CRMs are dispatched directly to the end-user in accordance with ISO requirements. Please check regulation status before ordering; additional fees can apply.

Product Overview

Dihydrocodeine (CRM) (Item No. ISO60142) is a certified reference material that is structurally categorized as an opioid. Dihydrocodeine is a weak μ -opioid receptor agonist ($K_i = 0.3 \mu\text{M}$) with analgesic, antitussive, and antidiarrheal applications. {32097,32099} It is metabolized by cytochrome P450 isoforms 2B6, 2D6, and 3A4, as well as uridine 5'-diphosphoglucuronosyltransferase, to products that are also active μ -opioid agonists. {32096,32098,32085} This product is intended for research and forensic applications.

FP/03/24

For research laboratory use only – Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

Contact Bertin Bioreagent

Parc d'Activit s du Pas du Lac
10 bis avenue Amp re
78180 Montigny le Bretonneux - France
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>
<https://www.bertin-bioreagent.com>