

Product Information



RXR $\hat{\pm}$ (Retinoid X receptor $\hat{\pm}$) Monoclonal Antibody

Cat No: G01034 - 100 $\hat{\mu}$ L

General Data

Shipping:	wet ice
Formulation:	Liquid Ascite, does not contain any preservative therefore avoid repeat freezing and thawing cycles
Host:	Mouse
Antigen:	recombinant protein, epitope corresponding to DE-region of mouseRXR α
Clone:	3A2
Isotype:	IgG1, κ



Application(s):	ELISA Western Blot Immunocytochemistry Immunohistochemistry Immunoprecipitation Gel shift Recommended dilutions: 1/500-1/5000
------------------------	---

Specificity: mouse and human RXR α . No cross reactivity on RXR β and RXR γ .

Scientific Literature

Rochette-Egly C, Lutz Y, Pfister V, Heyberger S, Scheuer I, Chambon P, Gaub MP. Detection of retinoid X receptors using specific monoclonal and polyclonal antibodies. *Biochem Biophys Res Commun.* 1994 Oct 28;204(2):525-36.

FP/04/24

For research laboratory use only – Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

Contact Bertin Bioreagent

Parc d'Activités du Pas du Lac
10 bis avenue Ampère
78180 Montigny le Bretonneux - France
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>
<https://www.bertin-bioreagent.com>