




Acetaminophen-d₄ (CRM)

Cat No: 22646 - 100 µg

General Data

Shipping: dry ice

Formulation: A 100 µg/ml solution in methanol



Special advice: Reference Materials RMs and Certified Reference Materials CRMs are dispatched directly to the end-user in accordance with ISO requirements. Please check regulation status before ordering; additional fees can apply.

Product Overview

Acetaminophen-d₄ (CRM) (Item No. 22646) is a certified reference material intended for use as an internal standard for the quantification of acetaminophen (Item Nos. 32839 22629) by GC- or LC-MS. Acetaminophen is categorized as an analgesic and antipyretic agent.^{42476,42477} Acetaminophen has been found as an adulterant in seized heroin samples.^{59755} Formulations containing acetaminophen have been associated with overdoses.^{54705} This product is intended for analytical forensic applications. This product is also available as a general research tool (Item No. 18245).

FP/26/24

For research laboratory use only – Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

Contact Bertin Bioreagent

Parc d'Activités du Pas du Lac
10 bis avenue Ampère
78180 Montigny le Bretonneux - France

Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>

<https://www.bertin-bioreagent.com>