



Anti-EP₂ Receptor:SureLight™ R-PE

Cat No: 16684 - 50 µg

General Data

Shipping: wet ice

Formulation: A lyophilized powder



Application(s): FC, FACS, Protein Microarrays

Product Overview

Anti-EP₂ (Prostaglandin E₂) Receptor conjugated to SureLight™ R-PE Antibody: Rabbit anti-EP₂ (Prostaglandin E₂) Receptor
Reactivity: Reactive with human, murine, and rat EP₂ receptors (other species not tested); Not reactive with EP₁, EP₃, and EP₄ receptors. Uses: Flow cytometry, cell-based assays, microarrays, and microplate applications

FP/22/24

For research laboratory use only – Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

Contact Bertin Bioreagent

Parc d'Activités du Pas du Lac
10 bis avenue Ampère
78180 Montigny le Bretonneux - France
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>
<https://www.bertin-bioreagent.com>