



Goat Anti-Rabbit IgG (H&L) Conjugated to SureLight[®] 488

Cat No: 27093 - 1 mg

General Data

Shipping: wet ice

Isotype: IgG



Application(s): Microscopy

Product Overview

Goat anti-Rabbit IgG (H&L) conjugated to SureLight[™] 488 reacts with the heavy and light chains of all rabbit immunoglobulins and subclasses. This antibody has no reactivity of the bovine, goat, human, mouse and rat immunoglobulins or serum proteins. Goat serum was obtained from healthy animals of US origin Excitation max. λ : 496 nm Emission max. λ : 517 nm Use: Microscopy

FP/04/24

For research laboratory use only – Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

Contact Bertin Bioreagent

Parc d'Activités du Pas du Lac
10 bis avenue Ampère
78180 Montigny le Bretonneux - France

Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>

<https://www.bertin-bioreagent.com>