




JNK1/2/3 (Phospho-Thr¹⁸³/Tyr¹⁸⁵) Rabbit Monoclonal Antibody

Cat No: 35906 - 100 µl

General Data

Shipping:	dry ice
Formulation:	100 µl of protein A-affinity purified monoclonal antibody
Host:	Rabbit
Antigen:	Peptide corresponding to human JNK1/2/3 (phospho-Thr183/Tyr185)
Clone:	RM464
Isotype:	IgG
Application(s):	IHC, WB



FP/29/24

For research laboratory use only – Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

Contact Bertin Bioreagent

Parc d'Activités du Pas du Lac
10 bis avenue Ampère
71810 Montigny le Bretonneux - France

Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>

<https://www.bertin-bioreagent.com>