

Product Information



5 α -Androstan-3 β ,17 β -diol Polyclonal Antibody

Cat No: C03155 - 500 dtn

General Data

Shipping: wet ice

Formulation: 1/2 glyceroled liquid

Host: Rabbit

Antigen: 5 α -Androstan-3 β ,17 β -diol-15-CMO-BSA

Application(s): RIA

Specificity: % CROSS REACTIVITY for B/Bo = 50%
C19
4 Androsten 17 -ol, 3-one (Testosterone) 3
5 Androstan 17 -ol, 3-one (Dihydrotestosterone : DHT) 5
5 Androstan 3 , 17 -diol 100
0,1
4 Androsten 3, 17-dione (4 Androstenedione) 0,6
5 Androsten 3 , 17 -diol (Andros



FP/22/24

For research laboratory use only – Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

Contact Bertin Bioreagent

Parc d'Activit s du Pas du Lac
10 bis avenue Amp re
78180 Montigny le Bretonneux - France
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>
<https://www.bertin-bioreagent.com>