

Product Information



6-keto Prostaglandin F1 α Polyclonal Antibody (P)

Cat No: C03051 - 500 dtn

General Data

Shipping:	wet ice
Formulation:	This vial contains 500 determinations of antibody in 1/2 glycerated liquid for a total volume of 25 μ L.
Host:	Rabbit
Antigen:	6KF1 α -BSA
Application(s):	RIA
Specificity:	% CROSS REACTIVITY for B/Bo = 50 % A1 0,1 B1 < 0,01 D1 0,4 E1 4,7 DH E1 19,9 K E1 < 0,06 DHK E1 < 0,01 F1 0,4 DHF1 1,3 K F1 < 0,01 DHKF1 < 0,01 6 keto E1 7 6 keto F1 100 6, 15 diketo F1 0,2 DH 6, 15 diketo F1 0,3 2, 3 dinor 6 keto F1 120



FP/22/24

For research laboratory use only – Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

Contact Bertin Bioreagent

Parc d'Activit s du Pas du Lac
10 bis avenue Amp re
78180 Montigny le Bretonneux - France
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>
<https://www.bertin-bioreagent.com>