

Product Information



Substance P Polyclonal Antibody IP

Cat No: C03082 - 500 dtn

General Data

Shipping: Dry Ice

Formulation: 1/2 glycerolated liquid

Host: Rabbit

Antigen: Substance P-OVA



Application(s): RIA

Specificity: % CROSS REACTIVITY for B/Bo = 50%

- Arg-Pro-Lys-Pro-Gln-Gln-Phe-Phe-Gly-Leu-Met-NH₂ (Subst P) 100
- Phe-Phe-Gly-Leu-Met-NH₂ 15
- Phe-Phe-Gly-Leu-Met-NH₂ 90
- Gln-Gln-Phe-Phe-Gly-Leu-Met-NH₂ 100
- Pro-Gln-Gln-Phe-Phe-Gly-Leu-Met-NH₂ 160
- Lys-Pro-Gln-Gln-Phe-Phe-Gly-L

FP/22/24

For research laboratory use only – Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

Contact Bertin Bioreagent

Parc d'Activités du Pas du Lac
10 bis avenue Ampère
78180 Montigny le Bretonneux - France
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>
<https://www.bertin-bioreagent.com>