

Product Information



TIF1 β Monoclonal Antibody

Cat No: G01100 - 100 μ L

General Data

Shipping:	wet ice
Formulation:	Liquid Ascite, does not contain any preservative therefore avoid repeat freezing and thawing cycles
Host:	Mouse
Antigen:	recombinant mouse TIF1 β as a fusion protein.
Clone:	1A9
Isotype:	IgG2b, κ



Application(s):	ELISA
	Western Blot
	Immunofluorescence
	Immunoprecipitation
	Recommended dilutions: 1/500-1/5000

Specificity:	human and mouse TIF1 β . Epitope localized between aa 123 and 834 of mouse TIF1 β . No cross reactivity on TIF1 α .
---------------------	--

Scientific Literature

CAMMAS F., OULAD ABDELGHANI M., VONESCH J.L., CHAMBON P. and LOSSON R. (2002). J. Cell. Sci. 115: 3439-3448

NIELSEN A.L., ORTIZ J.A., YOU J., OULAD ABDELGHANI M., KHECHUMIAN R., GANSMULLER A., CHAMBON P. AND LOSSON R. (1999). EMBO J. 18: 6385-6395

REMBOUTSIKA E, LUTZ Y, GANSMULLER A, VONESCH JL, LOSSON R, CHAMBON P. J. (1999). Cell. Sci. 112: 1671-1683

FP/03/24

For research laboratory use only – Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

Contact Bertin Bioreagent

Parc d'Activités du Pas du Lac
10 bis avenue Ampère
78180 Montigny le Bretonneux - France
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>
<https://www.bertin-bioreagent.com>