

Product Information



SAP130 \hat{A} Monoclonal Antibody

Cat No: G01062 - 100 \hat{A} μ L

General Data

Shipping:	wet ice
Formulation:	Liquid Ascite, does not contain any preservative therefore avoid repeat freezing and thawing cycles
Host:	Mouse
Antigen:	Synthetic peptide coupled to Ovalbumin
Clone:	1B10
Isotype:	IgG1, κ



Application(s):	ELISA Western Blot Immunofluorescence Recommended dilutions: 1/500 – 1/5000
------------------------	--

Specificity: human and mouse SAP130

Scientific Literature

Brand M, Moggs JG, Oulad-Abdelghani M, Lejeune F, Dilworth FJ, Stevenin J, Almouzni G, Tora L. UV-damaged DNA-binding protein in the TFIIIC complex links DNA damage recognition to nucleosome acetylation. EMBO J. 2001 Jun 15;20(12):3187-96

FP/03/24

For research laboratory use only \hat{A} Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

Contact Bertin Bioreagent

Parc d'Activit s du Pas du Lac
10 bis avenue Amp re
78180 Montigny le Bretonneux - France
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>
<https://www.bertin-bioreagent.com>