

Product Information



TAF10 \hat{A} Monoclonal Antibody

Cat No: G01076 - 100 \hat{A} μ L

General Data

Shipping:	wet ice
Formulation:	Liquid Ascite, does not contain any preservative therefore avoid repeat freezing and thawing cycles
Host:	Mouse
Antigen:	Full length human TAFI30 as a fusion protein
Clone:	2F4
Isotype:	IgG1, κ



Application(s):	ELISA Western Blot (endogenous and recombinant proteins) Immunocytofluorescence In immunoprecipitation works efficiently, but cannot be eluted with the peptide corresponding to the epitope Recommended dilutions: 1/500-1/5000
Specificity:	Human, does not work on mouse. Recognises the GGPGAGAAAGGTG epitope from aa. 50 to aa. 62

Scientific Literature

- Jacq, X., Brou, C., Lutz, Y., Davidson, I., Chambon, P. and Tora, L. (1994) Cell, 79, 107-117
- Wieczorek, E., Brand, M., Jacq, X. and Tora, L. (1998) Nature, 393, 187-191
- Metzger, D., Scheer, E., Soldatov, A., and Tora, L. (1999) EMBO J., 18, 4823-4834.

FP/03/24

For research laboratory use only \hat{A} Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

Contact Bertin Bioreagent

Parc d \hat{A} Activit \hat{A} s du Pas du Lac
10 bis avenue Amp \hat{A} re
78180 Montigny le Bretonneux - France
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>
<https://www.bertin-bioreagent.com>