

Product Information



TIF1 α Monoclonal Antibody


Cat No: G01099 - 100 μ L

General Data

Shipping:	wet ice
Formulation:	Liquid Ascite, does not contain any preservative therefore avoid repeat freezing and thawing cycles
Host:	Mouse
Antigen:	recombinant protein
Clone:	1E8
Isotype:	IgG1, κ

Application(s):	ELISA
	Western Blot
	Immunofluorescence
	Recommended dilutions: 1/500-1/5000

Specificity: human and mouse TIF1 α



Scientific Literature

REMBOUTSIKA E, LUTZ Y, GANSMULLER A, VONESCH JL, LOSSON R, CHAMBON P. J. (1999). Cell. Sci. 112: 1671-1683

NIELSEN A.L., ORTIZ J.A., YOU J., OULAD ABDELGHANI M., KHECHUMIAN R., GANSMULLER A., CHAMBON P. AND LOSSON R. (1999). EMBO J. 18: 6385-6395

FP/22/24

For research laboratory use only – Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

Contact Bertin Bioreagent

Parc d'Activit s du Pas du Lac
10 bis avenue Amp re
78180 Montigny le Bretonneux - France

Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>

<https://www.bertin-bioreagent.com>