

Product Information



TBPÂ Monoclonal Antibody

Cat No: G01093 - 100 ÂµL

General Data

Shipping:	wet ice
Formulation:	Liquid Ascite, does not contain any preservative therefore avoid repeat freezing and thawing cycles
Host:	Mouse
Antigen:	Full length recombinant human TBP.
Clone:	2C1
Isotype:	IgG2b, κ



Application(s):	ELISA Western Blot (endogenous and recombinant proteins) Immunocytofluorescence In immunoprecipitation works efficiently and the IP-ed complexes can be eluted with the GAMTPGIPIFSPMMPYG epitope peptide Recommended dilutions: 1/500-1/5000
------------------------	---

Specificity: Human, mouse, chicken

Scientific Literature

- Jacq, X., Brou, C., Lutz, Y., Davidson, I., Chambon, P. and Tora, L. (1994) Cell, 79, 107-117
Wieczorek, E., Brand, M., Jacq, X. and Tora, L. (1998) Nature, 393, 187-191.

FP/22/24

For research laboratory use only - Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

Contact Bertin Bioreagent

Parc d'Activit s du Pas du Lac
10 bis avenue Amp re
78180 Montigny le Bretonneux - France
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>
<https://www.bertin-bioreagent.com>