

Product Information



Heterochromatin Protein 1 α (HP-1 α) Monoclonal Antibody

Cat No: G01046 - 100 μ L

General Data

Shipping:	wet ice
Formulation:	Liquid Ascite, does not contain any preservative therefore avoid repeat freezing and thawing cycles
Host:	Mouse
Antigen:	Recombinant mouse HP1 α 0
Clone:	2G9
Isotype:	IgG1, κ



Application(s):	ELISA Western blot Immunoprecipitation ChIP For Immunohistochemistry and Immunofluorescence we recommend the clone 2HP-1H5 Recommended dilution : 1/500 – 1/ 5000
------------------------	--

Specificity: Hp1 α protein from mouse and human, epitope between aa 67 and 119

Scientific Literature

REMBOUTSIKA E, LUTZ Y, GANSMULLER A, VONESCH JL, LOSSON R, CHAMBON P. J. (1999). Cell.Sci. 112: 1671-1683

NIELSEN A.L., ORTIZ J.A., YOU J., OULAD ABDELGHANI M., KHECHUMIAN R., GANSMULLER A.,CHAMBON P. AND LOSSON R. (1999). EMBO J. 18: 6385-6395

FP/22/24

For research laboratory use only – Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

Contact Bertin Bioreagent

Parc d'Activités du Pas du Lac
10 bis avenue Ampère
78180 Montigny le Bretonneux - France
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>
<https://www.bertin-bioreagent.com>