

# Product Information



## PPAR $\beta$ 2 (Phospho-S112) Monoclonal Antibody

Cat No: G01025 - 100  $\mu$ L

### General Data

<b>Shipping:</b>	wet ice
<b>Formulation:</b>	Liquid Ascite, does not contain any preservative therefore avoid repeat freezing and thawing cycles
<b>Host:</b>	Mouse
<b>Antigen:</b>	Phosphorylated (serine 112) synthetic peptide coupled to ovalbumin: aa 107-121 of human PPAR $\beta$ 2, identical to PPAR $\beta$ 1 (KVEPA[phosphoS]PPYYSEKTC). No homology to PPAR $\alpha$ and PPAR $\beta$
<b>Clone:</b>	1F5
<b>Isotype:</b>	IgG1, $\kappa$
<b>Application(s):</b>	ELISA Western Blot Immunofluorescence recommended dilution: 1/500-1/5000
<b>Specificity:</b>	human and mouse phosphorylated PPAR $\beta$ 1, 2



FP/22/24

For research laboratory use only – Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

Contact Bertin Bioreagent

Parc d'Activit s du Pas du Lac  
10 bis avenue Amp re  
78180 Montigny le Bretonneux - France  
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>  
<https://www.bertin-bioreagent.com>