



## TFIIH p62 (TFIIH basal transcription factor p62) $\hat{A}$ Monoclonal Antibody

Cat No: G01098 - 100  $\hat{A}$  $\mu$ L

### General Data

|                     |   |
|---------------------|---|
| <b>Shipping:</b>    | wet ice   |
| <b>Formulation:</b> | Liquid Ascite, does not contain any preservative therefore avoid repeat freezing and thawing cycles |
| <b>Host:</b>        | Mouse   |
| <b>Antigen:</b>     | Recombinant human TFIIHp62.   |
| <b>Clone:</b>       | 3C9   |
| <b>Isotype:</b>     | IgG1, $\kappa$  |



|                        |  |
|------------------------|--|
| <b>Application(s):</b> | ELISA<br>Western Blot<br>Immunoprecipitation<br>Immunocytochemistry<br>Recommended dilutions: 1/500-1/5000 |
|------------------------|--|

**Specificity:** human TFIIHp62

### Scientific Literature

Coin F, Bergmann E, Tremeau-Bravard A, Egly JM. Mutations in XPB and XPD helicases found in xeroderma pigmentosum patients impair the transcription function of TFIIH. *EMBO J.* 1999 Mar 1;18(5):1357-66.

Coin F, Frit P, Viollet B, Salles B, Egly JM. TATA binding protein discriminates between different lesions on DNA, resulting in a transcription decrease. *Mol Cell Biol.* 1998 Jul;18(7):3907-14.

Rosignol M, Kolb-Cheyne I, Egly JM. Substrate specificity of the cdk-activating kinase (CAK) is altered upon association with TFIIH. *EMBO J.* 1997 Apr 1;16(7):1628-37.

FP/03/24

**For research laboratory use only  $\hat{A}$  Not for human diagnostic use.**

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

**Contact Bertin Bioreagent**

Parc d'Activit s du Pas du Lac  
10 bis avenue Amp re  
78180 Montigny le Bretonneux - France  
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>  
<https://www.bertin-bioreagent.com>