



## TAFII 68Å Monoclonal Antibody

Cat No: G01090 - 100 ÅµL

### General Data

<b>Shipping:</b>	wet ice
<b>Formulation:</b>	Liquid Ascite, does not contain any preservative therefore avoid repeat freezing and thawing cycles
<b>Host:</b>	Mouse
<b>Antigen:</b>	The RNA-binding region (RNP) of human TAFII68 as a GST-fusion protein
<b>Clone:</b>	2B10
<b>Isotype:</b>	IgG1, Å



<b>Application(s):</b>	ELISA Western Blot (endogenous and recombinant proteins) Immunocytofluorescence Immunoprecipitation Recommended dilutions: 1/500-1/5000
------------------------	---

**Specificity:** Human, other species have not been tested.

### Scientific Literature

Bertolotti, A., Lutz, Y., Chambon, P. and Tora, L. (1996) EMBO J. 15, 5022-5031.

Bertolotti A., Melot T., Acker J., Vigneron M., Delattre O. and Tora L. (1998) Mol. Cell. Biol. 18, 1489-1497.

FP/04/24

**For research laboratory use only – Not for human diagnostic use.**

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

**Contact Bertin Bioreagent**

Parc d'Activit s du Pas du Lac  
10 bis avenue Amp re  
78180 Montigny le Bretonneux - France  
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>  
<https://www.bertin-bioreagent.com>