

# Product Information



## TAF5<sup>Δ</sup> Monoclonal Antibody

Cat No: G01085 - 100  $\mu$ L

### General Data

<b>Shipping:</b>	wet ice
<b>Formulation:</b>	Liquid Ascite, does not contain any preservative therefore avoid repeat freezing and thawing cycles
<b>Host:</b>	Mouse
<b>Antigen:</b>	Synthetic peptide coupled to ovalbumin
<b>Clone:</b>	2D2
<b>Isotype:</b>	IgG1, $\kappa$



<b>Application(s):</b>	ELISA Western Blot (endogenous and recombinant proteins) Immunocytofluorescence Immunoprecipitation Recommended dilutions: 1/500-1/5000
------------------------	---

**Specificity:** Human, does not work on mouse.

### Scientific Literature

Dubrovskaya, V., Lavigne, A.-C., Davidson, I., Acker, J., Staub, A. and Tora L. (1996). EMBO J. 15, 3702-3712.

Wieczorek, E., Brand, M., Jacq, X. and Tora, L. (1998) Nature, 393, 187-191.

FP/04/24

**For research laboratory use only – Not for human diagnostic use.**

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

**Contact Bertin Bioreagent**

Parc d'Activit s du Pas du Lac  
10 bis avenue Amp re  
78180 Montigny le Bretonneux - France  
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>  
<https://www.bertin-bioreagent.com>