



## TRRAP Monoclonal Antibody

Cat No: G01043 - 100 µL

### General Data

<b>Shipping:</b>	wet ice
<b>Formulation:</b>	Liquid Ascite, does not contain any preservative therefore avoid repeat freezing and thawing cycles
<b>Host:</b>	Mouse
<b>Antigen:</b>	Synthetic peptide coupled to Ovalbumin
<b>Clone:</b>	2D5
<b>Isotype:</b>	IgG1, κ



<b>Application(s):</b>	ELISA Western Blot: natural and recombinant protein Immunofluorescence Immunoprecipitation Chromatin IP Recommended dilutions: 1/500 – 1/5000
------------------------	--

**Specificity:** Human , mouse

### Scientific Literature

Helmlinger D, Hardy S, Sasorith S, Klein F, Robert F, Weber C, Miguët L, Potier N, Van-Dorsseleer A, Wurtz JM, Mandel JL, Tora L, Devys D. Ataxin-7 is a subunit of GCN5 histone acetyltransferase-containing complexes. *Hum Mol Genet.* 2004 Jun 15;13(12):1257-65. Epub 2004 Apr 28

Robert F, Hardy S, Nagy Z, Baldeyron C, Murr R, Dery U, Masson JY, Papadopoulou D, Herceg Z, Tora L. The Transcriptional Histone Acetyltransferase Cofactor TRRAP Associates with the MRN Repair Complex and Plays a Role in DNA Double-Strand Break Repair. *Mol Cell Biol.* 2006 Jan;26(2):402-12.

FP/04/24

**For research laboratory use only – Not for human diagnostic use.**

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

**Contact Bertin Bioreagent**

Parc d'Activités du Pas du Lac  
10 bis avenue Ampère  
78180 Montigny le Bretonneux - France  
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>  
<https://www.bertin-bioreagent.com>