



## Prostaglandin F2 $\alpha$ Polyclonal Antibody (P)

Cat No: C03054 - 500 dtn

### General Data

**Shipping:** wet ice

**Formulation:** This vial contains 500 determinations of antibody in 1/2 glycerated liquid for a total volume of 33.33  $\mu$ L.

**Host:** Rabbit

**Antigen:** PGF2 $\alpha$ -BSA

**Application(s):** RIA

**Specificity:** % CROSS REACTIVITY for B/Bo = 50 %

A1 <0.01  
B1 <0.01  
D1 0.4  
E1 0.03  
KE1 <0.01  
F1 12  
DHF1 0.4  
K F1 <0.01  
DHKF1 <0.01  
6-keto F1 0.04  
A2 <0.01  
D2 3  
E2 0.03  
K E2 <0.01  
TXB2 <0.01  
F2 100  
DHF2 3  
K F2 <0.01  
DHK F2 <0.01  
19 OH F2 2  
A2 <0.01  
D2 3  
E2 0.03  
K E2 <0.0



FP/22/24

For research laboratory use only – Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

Contact Bertin Bioreagent

Parc d'Activit s du Pas du Lac  
10 bis avenue Amp re  
78180 Montigny le Bretonneux - France  
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>  
<https://www.bertin-bioreagent.com>