

# Product Information



## Ultra pure water ELISA Grade

Cat No: A07001 - 1L

### General Data

**Shipping:** room temp.



**Application(s):** Preparation of all EIA buffers and reagents.

### Product Overview

The Ultra pure water is deionised demineralised water. It was produced using a reverse osmosis, ion exchange and organic scavenging to ensure the highest level of purity. To be used to prepare all EIA buffers and reagent.

FP/22/24

**For research laboratory use only – Not for human diagnostic use.**

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

**Contact Bertin Bioreagent**

Parc d'Activités du Pas du Lac  
10 bis avenue Ampère  
78180 Montigny le Bretonneux - France

Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>

<https://www.bertin-bioreagent.com>