

Product Information



Prostaglandin E2 Monoclonal Antibody

Cat No: C03011 - 200 µg

General Data

Shipping:	Dry Ice
Formulation:	This vial contains 200 µg of antibody in 1/2 glycerolated liquid
Host:	Mouse blb C myelome SP 2,0 Ag
Antigen:	PGE2α-Thyroglobulin
Clone:	8,00E+57

Application(s): EIA, RIA

Specificity:	% CROSS REACTIVITY for B/Bo = 50 %
	A1 <0.06
	D1 0.03
	DH E1 0.1
	DHK E1 <0.01
	K F1 <0.01
	6-keto E1 10
	8-iso PGF2 0.63
	B2 0.01
	E2 100
	K E2 0.08
	F2 1.5
	DH F2 <0.01
	DHK F2 <0.01
	19 OH E2 0.05
	8-iso E2 8
	15-epi E2 0.06
	B1 <0.01
	E1 70
	K E1 0.07
	F1 DHK F1 <0.01
	6-



FP/03/24

For research laboratory use only – Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

Contact Bertin Bioreagent

Parc d'Activités du Pas du Lac
10 bis avenue Ampère
78180 Montigny le Bretonneux - France
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>
<https://www.bertin-bioreagent.com>