

Product Information



6-keto Prostaglandin E1 Polyclonal Antibody

Cat No: C03019 - 500 dtn

General Data

Shipping: wet ice

Formulation: This vial contains 500 determinations of antibody in 1/2 glycerated liquid for a total volume of 1,25 µL.

Host: Rabbit

Antigen: 6-keto PGE1-BSA

Application(s): R/A

Specificity: % CROSS REACTIVITY for B/Bo = 50 %
A1 0,01%
B1 0,4%
D 0,15%
E1 0,08%
DH E1 0,03%
K E1 < 0,01%
DHK E1 < 0,01%
F1 0,01%
DHF1 < 0,01%
19 OH E1 methyl ester < 0,01%
6 keto E1 100%
6 keto F1 0,07%
6, 15 diketo F1 0,01%
DH 6, 15 diketo F1 0,16%
A2 0,26%
B2 1%
D



FP/03/24

For research laboratory use only – Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

Contact Bertin Bioreagent

Parc d'Activités du Pas du Lac
10 bis avenue Ampère
78180 Montigny le Bretonneux - France
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>
<https://www.bertin-bioreagent.com>