



## Mouse Anti-Fluorescein IgG:SureLight® R-PE

Cat No: 16691 - 100 µg

### General Data

<b>Shipping:</b>	wet ice
<b>Formulation:</b>	Lyophilized powder
<b>Clone:</b>	027-11741
<b>Isotype:</b>	IgG2A



**Application(s):** FC, FACS, Protein Microarrays

### Product Overview

Mouse anti-Fluorescein IgG2a (clone 027-11741) conjugated to SureLight® R-PE Antibody: Mouse anti-Fluorescein (FITC) IgG2a (clone 027-11741) Dye: SureLight® R-phycoerythrin Excitation max.  $\lambda$ : 565-498 nm Emission max.  $\lambda$ : 578 nm Uses: Flow cytometry, cell-based assays, microarrays, and microplate applications Clone Number: 027-11741 Isotype: IgG2A

FP/22/24

**For research laboratory use only – Not for human diagnostic use.**

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

**Contact Bertin Bioreagent**

Parc d'Activités du Pas du Lac  
10 bis avenue Ampère  
78180 Montigny le Bretonneux - France  
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>  
<https://www.bertin-bioreagent.com>