

Product Information




Cocaine-d₃ (CRM)

Cat No: 19423 - 1 mg

General Data

Shipping: dry ice

Formulation: A 100 µg/ml or 1 mg/ml solution in acetonitrile



Special advice: Reference Materials RMs and Certified Reference Materials CRMs are dispatched directly to the end-user in accordance with ISO requirements. Please check regulation status before ordering; additional fees can apply.

Product Overview

Cocaine-d₃ (CRM) (Item No. 19423) is intended for use as an internal standard for the quantification of cocaine (Item Nos. ISO60176

16186) by GC- or LC-MS. Cocaine is an analytical reference standard that is structurally categorized as a tropane. It modulates monoamine uptake and brain reward pathways with a high risk of dependence. Cocaine is regulated as a Schedule II drug in the United States. This product is intended for research and forensic applications.

FP/03/24

For research laboratory use only – Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

Contact Bertin Bioreagent

Parc d'Activités du Pas du Lac
10 bis avenue Ampère
78180 Montigny le Bretonneux - France

Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>

<https://www.bertin-bioreagent.com>